



# MarketLINK Scientific

## FastPette Pipette Controller



### Labnet FastPette Pipette Controller

**Complete Control**—aspiration speed can be set to high or low for fast filling of larger pipettes or precision measurement of small volumes.

**Protection and Safety** - a unique safety valve prevents overfilling and a PTFE filter protects samples from contamination. The filter, nosepiece and pipette holder are easily removed for autoclaving.

**Ready for Action**—when fully charged the NiMH battery enables up to 8 hours of continuous use. A low battery light indicates when the unit requires charging.



Catalogue # P2000

**\$279.00**

The FastPette is supplied with PTFE filter, power supply/charger, bench stand and wall hanger.



Contact us to find out more about the Discovery 4 Pack Starter Kit and the complete line of Labnet pipettes and supplies!

MarketLINK Scientific  
919 Fraser Drive Unit 11, Burlington, ON L7L 4X8  
P: 866-336-3313  
F: 289-288-0122

Applicable taxes not included. Prices subject to change.

For more information contact: Allison Bureau  
allison@marketlinkscientific.com  
P: 866-336-3313x3364  
www.marketlinkscientific.com